



Genomic Unity™ Test

Detecting tandem repeat expansions using WGS

Genomic Unity™ test detects and reports tandem repeat expansions in selected, characterized regions with known pathogenic associations.

What role do tandem repeat expansions play in genetic disease?


There are dozens of tandem repeat expansions that are known to be pathogenic. Some, like the CTG repeat expansion in the DMPK gene, lead to predominantly adult-onset neurological disorders like myotonic dystrophy type I. Others, like the CGG repeat expansion within the 5' UTR of the FMR1 gene which causes Fragile X syndrome, affect the individual early in life and are routinely tested for as part of the rare disease diagnostic odyssey.

At each pathogenic locus, alleles are classified as normal/intermediate, premutation or full mutation based on the number of repeats, with the ranges varying for each locus. For example, in the case of FMR1, alleles up to 44 repeats in length are considered normal. Alleles between 55 and 200 repeats in length are considered premutation and are at a high risk of having the number of repeats further expand into the full mutation range of >200 repeats when transmitted from a parent to child. Until recently, detecting repeat expansions has required the use of PCR or southern blot analysis, usually employed to interrogate a single targeted gene. Panel and exome-based tests will miss these genetic changes.

How does Variantyx detect tandem repeats?

Variantyx uses whole genome sequencing (WGS) technology to provide consistent, comprehensive coverage of the entire genome. To analyze tandem repeats in known pathogenic loci, three separate paired-end read strategies are used, all within a single assay.

The first strategy focuses on short range variants that are less than 50 repeats in length. Here de novo assembly of spanning reads is used as the full length of the repeat is contained within either the R1 or R2 read, with uniquely mappable flanking sequences.



The second strategy focuses on medium range variants that span from 50 to 200 repeats in length, with the upper limit determined by the sequencing insert size of 550 bp. Here alignment and counting of repeats in anchored reads is used. With

anchored reads, one member of the pair, either the R1 or R2 read, contains only repeat sequence while the other member contains partial repeat sequence and partial uniquely mappable sequence.

The final strategy focuses on long range variants that are >200 repeats in length. Here counting and statistical normalization of reads containing only repeat sequence is used to estimate the repeat length. Combining the three different methods, repeat length is calculated with good specificity up 200 repeats. Alleles with repeat lengths exceeding this threshold represent high-confidence estimates that should be independently confirmed by an orthogonal technology.

What tandem repeats are analyzed?

The following table shows analysis results from a 50 year old male patient with suspected myotonic dystrophy type II. Testing confirmed the myotonic dystrophy type II diagnosis.

Gene Repeat: Allele count	Disorder Allele interpretation	Gene Repeat: Allele count	Disorder Allele interpretation
AFF2 CCG: 12, -	Fragile XE syndrome Normal	CACNA1A CAG: 11, 12	Spinocerebellar ataxia Normal
AFF3 CGG: 8, 8	FRA2A fragile site Normal	CNBP CCTG: 110, 127	Myotonic dystrophy type II Full mutation
AR CAG: 23, -	Spinal and bulbar muscular atrophy Normal	CSTB CCCCGCCCGCG: 3, 2	Myoclonus epilepsy Normal
ATN1 CAG: 19, 12	Dentatorubral-pallidoluysian atrophy Normal	DIP2B CGG: 12, 12	FRA12A fragile site Normal
ATXN1 CAG: 31, 32	Spinocerebellar ataxia Normal	DMPK CTG: 5, 11	Myotonic dystrophy type I Normal
ATXN10 ATTCT: 12, 15	Spinocerebellar ataxia Normal	FMR1 CGG: 30, -	Fragile X syndrome Normal
ATXN2 CAG: 22, 23	Spinocerebellar ataxia Normal	FXN GAA: 9, 19	Friedreich's ataxia Normal
ATXN3 CAG: 20, 24	Spinocerebellar ataxia Normal	HTT CAG: 17, 20	Huntington disease Normal
ATXN7 CAG: 10, 10	Spinocerebellar ataxia Normal	JPH3 CTG: 14, 16	Huntington disease-like 2 syndrome Normal
ATXN80S CTG: 15, 16	Spinocerebellar ataxia Normal	PPP2R2B CAG: 10, 17	Spinocerebellar ataxia Normal
C9ORF72 GGGGC: 2, 5	FTDALS1 Normal		

What are the benefits of Genomic Unity™ tandem repeat analysis?

With Genomic Unity™ there is no need for a separate sample and separate assay.

Because WGS provides comprehensive coverage of the entire genome, all sequence data necessary for detection of tandem repeat expansions is present. Variantyx's custom-built, validated computer algorithms analyze the data, identifying repeat expansions for more than 20 different loci.